



University of Essex

University of Essex

School of Electronic Engineering and Computer Science
University of Essex

Haptic-Controlled Robotic Arm Interface with Dobot Magician

The School of Computer Science and Electronic Engineering, University of Essex, CO4 3SQ, Essex, UK.



Prerequisites and Setup

- Required Software:
 - Visual Studio (with Python & C++ extensions)
 - Python 3.13
 - DobotStudio
 - Haptic Device Software
- Project Folder:
 - Includes source code, DLLs, .bat file
- Environment Setup:
 - Add project folder path to environment variables

System Overview

- Core Files:

- Master-Slave_Device.cpp: Haptic control, force feedback, pipe output


The “Master-Slave_Device.cpp” file manages the connection and control of the haptic device and sends data to the Python script controlling the Dobot Magician (DobotControl.py).

- DobotControl.py: Reads pipe data, controls Dobot

“DobotControl.py” handles data processing from the haptic device and controls the Dobot.

- Data Flow:

Haptic Device → C++ → Pipe → Python → Dobot




Running the Program

1. Launch the .bat file

After the set up is complete, you can run the program by clicking the .bat-file in the project folder.

2. Two terminal windows will open:

- **Haptic communication**
- **Dobot control**



Adjusting Dobot Movements

- Python Control Command:
- `dType.SetJOGCmd(api, 0, cmd, 0)`

- **Command Selection:** Determines the appropriate command (`cmd`) for the robot based on the processed input values and button states.
 - `cmd = 1, 2` : Move along the X-axis (positive/negative).
 - `cmd = 3, 4` : Move along the Y-axis (positive/negative).
 - `cmd = 5, 6` : Move along the Z-axis (positive/negative).
 - `cmd = 7, 8` : Rotate the catheter (clockwise/counterclockwise).
 - `valueb = 1, 4` : Drive the catheter in or out using digital I/O pins.
 - `valueb = 5, 10` : Lock/unlock the robot's movement to allow for catheter manipulation without moving the arm.

For the Magnetic Guidewire's Application Test:

- 'cmd' controls direction
- Modified for X-axis only
- Movement scaled down for precision (For the application test is 10)

```
# set coordinates for the dobot reading the values from the pipe.
if lock == 0:
    if valuex>0:
        cmd = 1
    elif valuex<=0:
        cmd = 2

    elif valuey>0:
        cmd = 3
    elif valuey<=0:
        cmd = 4
    elif valuez>0:
        cmd = 5
    elif valuez<=0:
        cmd = 6
    elif valueb == 8:
        cmd = 7
    elif valueb == 2:
        cmd = 8
    else:
        cmd = 0

#map the velocity of the arm according to the radius of the haptic devices input
vel = math.sqrt(pow(valuex,2)+pow(valuey,2)+pow(valuez,2))

# send the command to move the dobot
dType.SetJOGCmd(api, 0, cmd, 0) #Uncomment this for Lorans project
```

This is the part of the code that controls the Dobot movement.

Guidewire Manipulation

- Use buttons to move catheter in/out
- Direction may be reversed depending on USB port
- Press two buttons simultaneously to lock/unlock Dobot movement




Troubleshooting - Dobot

- **Issue:** Blue light, unresponsive Dobot

- **Fix:**

1. Open DobotStudio
2. Run Blockly program (EIO10 or EIO13)
3. Disconnect, close terminals
4. Relaunch the program



Troubleshooting - Haptic Device

- **Issue:** Device not detected
- **Fix:**
 1. Open 'Teneo' test tool
 2. Confirm detection
 3. Try different USB port
 4. Reboot or reinstall haptic software if needed